


PALM User Protocols

Immunofluorescence on frozen sections



PALM Laboratories
A unique service of Carl Zeiss



We make it visible.

Immunofluorescence



Some helpful tips before starting:

To prepare sections for non-contact laser capture microdissection (LCM) we generally recommend **PALM MembraneSlides**. If weak fluorescence or small objects must be detected the MembraneSlide 1.0 PET (Order No. 415101-4401-050) is recommendable. (PEN membrane on PALM MembraneSlides can lead to more background in some fluorescent filters due to its structure and stronger autofluorescent effects)



Note: To allow subsequent cutting and lifting a **coverslip** and standard mounting medium **must not be applied!** Freezing media like OCT or similar may be used for sectioning but should be kept to a minimum and have to be removed before laser cutting.

For collecting microdissected samples we recommend the special **PALM AdhesiveCaps**: AdhesiveCap 500 opaque (Order No. 415101-4400-250) or AdhesiveCap 500 clear (Order No. 415101-4400-255).



Tongue epithelium
(Rat, Anti-Cytokeratin-FITC, DAPI)

Immunofluorescence

Microtome cutting and antibody staining

- Tissue (e.g., frozen rat tongue) is cut as usual (8-10 μm).
- The frozen section is transferred from the blade to the warmer slide by cautious touching.
- **Dry** in the cryostat for **1 minute** at about **-20°C**.

Note: *Few minutes of longer drying or previous poly-L-lysine coating of the slide may improve the adhesion of sections during the following staining steps.*

- Subsequently **fix** and dehydrate the section in **ice-cold** pure **acetone** for **30 seconds**.
- Finally **air-dry** at room temperature for **2-3 minutes**. Keep refrigerated till further use.
- To reduce the necessary incubation volumes draw a hydrophobic line around the sections with a special pen (e.g., Dako Pen #S2002).
- **Rehydration** and **blocking** of unspecific binding is done by covering the section with a drop of ready-to-use Protein Block Serum-Free (DakoCytomation #X0909) for **15 minutes** at room temperature.
- Pour off protein block and dip once into **PBS** for short **washing**.
- Remove excess liquid from the slide by tapping on an absorbant surface.
- Incubate with **Ab/DAPI-solution** for **1 hour** at room temperature in a dark wet chamber.
- Pour off Ab/DAPI-solution and dip once into **PBS** for short **washing**.
- Remove excess liquid from the slide by tapping on an absorbant surface.
- Sections can now be viewed on the microscope and used for LCM. To stabilize the fluorescent signals a drop of PBS or VECTASHIELD may be added as cover but microdissection will only be possible without too much liquid.

Note: *FITC-fluorescence is very sensible to bleaching without protection (often not more than 30 seconds) and therefore fluorescent illumination time should be kept to a minimum. The „Freeze Mode“-function of the PALM RoboSoftware will be very helpful in this context.*

Preparation of Ab/DAPI-staining solution

Dilute labelled **antibody** (Ab) solution in **PBS** containing **0.5 $\mu\text{g/ml}$ DAPI** shortly before use (e.g., mouse monoclonal PCK-26, panCytokeratin-FITC; #ab112114, abcam; e.g., 1:25 - 1:100). The optimal dilution has to be found empirically for any tissue. Keep this solution in the dark till use and during the incubation.